



**FOR PEOPLE
FOR INNOVATION
FOR HEALTHCARE**

REINFORCED ALGINATE DRESSINGS WITH SILVER

Antibacterial. Natural alginate fibres reinforced with cellulose fibres. Soft, conformable and easy to remove.

KEY FEATURES:

- High absorbency
- Reinforced
- Fast effective antibacterial activity¹
- Intact removal due to enhanced wet strength²
- Enhanced conformability
- Moist wound environment promoting healing

INDICATIONS:

- Moderately and highly exuding wounds
- Superficial or deep wounds

Such as:

- Superficial or partial thickness burns
- Cavity wounds
- Venous lower leg ulcers
- Diabetic ulcers
- Post-operative wounds
- Pressure ulcers


This dressing is effective against: Gram positive and negative bacteria, such as E. Coli and MRSA

AVAILABLE FORMAT:

Roll form, cut dressings, packaged and sterilised.

Sizes³:

- | | |
|-------------|-----------|
| • 5x5cm | • 10x20cm |
| • 10x10cm | • 15x15cm |
| • 10x12.5cm | • 2x45cm |



¹ PAgTr-49

² PLUPVR-Z, PLUPVR-4, PLUAgPVP-1, PLUAgPVP-2

³ Other sizes are available upon request

+44 (0)2476 708 200 | enquiries@sfm-limited.com | www.sfm-limited.com

SFM Limited, Galaxy House, 31 Herald Way, Binley Industrial Estate, Coventry, CV3 2RQ, United Kingdom